## Jen-Unibond LV

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## Jen-Unibond EN

Jen-Unibond is a universal self-priming dental adhesive agent designed to bond the composites materials to enamel and dentin.

### •• ADVANTAGES

- Jen-Unibond possesses the very high bond values of any single-component resin bonding system. Jen-Unibond bond values to etched dentin after 24 hours are in the high 30 – 50 MPa range depending on the protocol of individual laboratories.
- Contains the natural resins providing an elastic quality to all composition, which creates a hybrid layer; that is able to flex. A flexible hybrid layer reduces the risk of microleakage by creating the layer that; is not rigid and resistant to fracture; the natural resins contain also the natural surfactants, which, allow for easier penetration of resins into the dentin tubules.
- The special co-polymerized solvent in combination with very small amount of ethyl alcohol carrier; facilitates the problem of quick evaporation of unwanted solvent, allowing simultaneously for excellent penetration of the resin in the dentinal tubule structure.

### •• INDICATIONS

- A bonding agent for composite (including flowable) restorations;
- Adhesive cavity varnish under amalgam restorations;
- Protective varnish for hypersensitive tooth areas.

### •• CONTRAINDICATIONS

- Use for patients with allergy to methacrylates;
- Direct or indirect pulp capping.

### •• INSTRUCTION FOR USE

#### ISOLATION

- Rubber dam is the recommended method of isolation.

#### ETCHING

- Apply etching gel by brushing it into the dentin. Extend the application on the enamel margins. Leave on for 15 – 20 seconds.
- If a relatively large area of the enamel is involved in the restoration, and in clinical situations where achieving maximum bonding strength to enamel is critical to good retention, the enamel should be thoroughly treated with the etching gel for 20 seconds.
- Rinse well, then dry lightly with a stream of air. Do not desiccate the dentin structure.

#### DENTIN-ENAMEL PRIMER/BONDING AGENT

- Apply primer /bonding agent with a brush in a manner similar to the application of the cavity varnish. Make sure the dentin is well wetted. Just placing primer over the dentin may not be enough. Apply one or two consecutive coats of Jen-Unibond to etched enamel and dentin. Do not rinse off. Keep the area dry; avoid saliva contamination. Light cure for 10 seconds.

- Follow with the application of the light cured restorative material or cement and cure thoroughly. The required curing time may depend on the intensity of the curing light, the thickness of the layer of restorative material or cement, its color and reactivity. Usually a curing time of 40 seconds is enough. Refer to manufacturer’s instruction for using the restorative materials.